

# Spatial Immune Profiling Assay Kit

14-Plex for Human FFPE Tissue and ZellScannerONE™

AbKT-3001-10RXN

ChipCytometry™ Assay Kits contain ready-to-use, reliable reagents and optimized protocols to enable researchers to obtain quick, robust data with the ChipCytometry platform. Each ChipCytometry Assay Kit contains pre-validated fluorescent antibodies and buffers for staining 10 samples. Assay Kit optimization involves extensive reagent testing on positive and negative controls and review by a board-certified pathologist. ChipCytometry Assay Kits can be customized by adding antibodies for proteins of interest after the final staining cycle.



**Ready-to-use:** Get your results quickly with this optimized assay kit. Assay kits include pre-validated antibodies with optimized dilutions and staining plan to accelerate your experiment.



**Reliable:** Reagents are selected from a variety of vendors and are highly validated for use in ChipCytometry assays. Each assay kit has been tested in multiple donors to validate reproducibility.



**Customizable:** Each assay kit can be customized to suit researchers needs. Expand upon the kit with pre-validated antibodies from our catalog to customize your panel.

## Ready-to-Use Kit

Spatial Immune Profiling Assay Kits contain ready-to-use, reliable reagents and optimized protocols for quick, robust data collection. Identify immune and tumor cells in human FFPE tissue.

Each kit contains 14 pre-validated fluorescent antibodies, nuclear counterstains, and buffers for staining 10 samples. Additional materials, including Chip Kits and Wash Buffer, are required for experiments using the ZellScannerONE and must be purchased separately.

- **Product Number:** AbKT-3001-10RXN
- **Sample Type:** Human FFPE tissue
- **Size:** 10 reactions



The Spatial Immune Profiling Assay Kit contains materials for 10 reactions. The kit includes pre-validated fluorescent antibodies, nuclear counterstains, antibody diluent, and staining protocol for human FFPE tissue.

Staining Protocol					
Cycle	Target	Filter Set	Antibody Volume	Diluent Volume	Incubation Time
1	Ki67	FS488	50 $\mu$ l	450 $\mu$ l	1 hr
	CD27	FS560	50 $\mu$ l		
	CD3	FSPerCP	50 $\mu$ l		
2	CD4	FS488	50 $\mu$ l	500 $\mu$ l	1 hr
	CD45	FSPerCP	50 $\mu$ l		
3	CD68	FS488	50 $\mu$ l	450 $\mu$ l	1 hr
	FoxP3	FS560	50 $\mu$ l		
	CD45RA	FSPerCP	50 $\mu$ l		
4	CD8	FS488	50 $\mu$ l	450 $\mu$ l	1 hr
	PD-1	FS560	50 $\mu$ l		
	EpCAM	FSPerCP	50 $\mu$ l		
5	Granzyme B	FS488	50 $\mu$ l	500 $\mu$ l	1 hr
	PD-L1	FS560	50 $\mu$ l		

The staining protocol for the Spatial Immune Profiling Assay Kit is accomplished in just 5 cycles. To customize your panel, add additional cycles selecting from our catalog of pre-validated antibodies or supplement with antibodies from your own inventory.

## Customizable Protocol

Get your results quickly with this optimized 14-plex protocol. Assay kits include pre-validated antibodies with optimized dilutions and staining plan to accelerate your experiment. This pre-optimized protocol enables researchers to get to results faster.

Each assay kit can be customized to suit individual project needs. Expand upon the kit with pre-validated antibodies from our catalog to customize your panel.



## Reliable Coverage

Spatially profile key immune cell phenotypes with the Spatial Immune Profiling Assay Kit. Antibodies in the kit are selected from a variety of vendors and are highly validated for use in ChipCytometry assays. Assay Kit optimization involves extensive testing on control FFPE tissues and staining review by a board-certified pathologist.

Identify cell phenotypes critical for your spatial profiling assay. Phenotypes include key T cell subtypes, macrophages, and epithelial cells. Characterize the cells in your FFPE tissues using the gating strategy below or define your own gating strategy to identify cell types of interest.

Markers											
Phenotype	CD45	EpCAM	CD3	CD4	CD8	CD27	CD68	FoxP3	Ki-67	PD-1	CD45RA
Leukocytes	✓										
Epithelial cells		✓									
T cells	✓		✓								
Macrophages	✓						✓				
T helper cells	✓		✓	✓							
T regulatory	✓		✓	✓				✓			
Proliferating T helper cells	✓		✓	✓					✓		
Exhausted T helper cells	✓		✓	✓						✓	
Central memory T helper cells	✓		✓	✓		✓					
Naive T helper cells	✓		✓	✓		✓					✓
Effector T helper cells	✓		✓	✓							✓
Effector memory T helper cells	✓		✓	✓							✓
T cytotoxic cells	✓		✓		✓						
Proliferating T cytotoxic cells	✓		✓		✓				✓		✓
Exhausted T cytotoxic cells	✓		✓		✓					✓	
Naive T cytotoxic cells	✓		✓		✓	✓					✓
Effector T cytotoxic cells	✓		✓		✓						✓
Effector memory CD8+ T cells	✓		✓		✓						

The Spatial Immune Profiling Kit enables spatial phenotyping of key immune populations and epithelial cells, including those listed in the table above. This is partial list of the 8,192 distinct phenotypes that can be identified with this panel. PD-L1 and Granzyme B phenotypes are not represented in this table. Additional phenotypes can be identified for different degrees of expression of single markers.



---

To request a quote, visit [canopybiosciences.com/chipcytometry](https://canopybiosciences.com/chipcytometry)

---

Canopy Biosciences® - A Bruker Company  
4340 Duncan Ave., Ste 220  
St. Louis, MO 63110

**For Research Use Only.** Not for use in diagnostic procedures.

© 2022 Canopy Biosciences – A Bruker Company. All rights reserved. Canopy Biosciences and the Canopy Biosciences logo are trademarks or registered trademarks of Canopy Biosciences – A Bruker Company, in the United States and/or other countries.

AUGUST 2022 V2