

#309 Seeing more of the brain with high-plex multi-omics imaging at single-cell spatial resolution in a large contiguous tissue section

Dwayne Dunaway¹, Zachary R. Lewis¹, Victoria Rachleff², Tushar Rane¹, Tiên Phan-Everson¹, Alyssa Rosenbloom¹, Shilah A. Bonnett¹, Giang Ong¹, Lidan Wu¹, Ashley Heck¹, Kimberly Young¹, Rustem Khafizov², David Ruff², Erin Piazza¹, Patrick Danaher¹, Aster Wardhani¹, Mithra Korukonda¹, Carl Brown¹, Christine Kang¹, John Lyssand¹, Jim Forrester¹, Trieu My Van¹, Vik Devgan¹, Jaemyeong Jung¹, Gary Geiss¹, Joseph M. Beechem¹, and C. Dirk Keene²
¹NanoString® Technologies, Seattle, WA, USA ²University of Washington, Department of Laboratory Medicine and Pathology, Seattle, WA, USA

Abstract


Neurophysiological and neurodegenerative processes occur over multiple spatial scales in the brain. While fine-scaled spatial molecular analyses of human brain function have vastly advanced in the last decade, these datasets often lack context of how microscopic features fit within the larger picture of macroscopic neuroanatomy. In the brain, critical cell functions and cell-to-cell communications take place over relatively large interstitial distances and far from the cell somas. The ability to explore protein and RNA-driven activities at high resolution, within the spatial context across a large environment and passing through multiple microenvironments, allows a comprehensive picture of brain biology to emerge and is applicable to open questions regarding brain development, activity, aging, and neurodegeneration.

Here we demonstrate a novel Larger Surface Area flow cell with > 1,600 mm² of imageable area. Fluidic and imaging conditions of Spatial Molecular Imager (SMI) are optimized to deliver consistent performance. Combining the Larger Surface Area flow-cell imaging capabilities with the human-specific > 68-plex CosMx™ SMI Human Neuroscience protein panel (Neural Cell Profiling and Alzheimer's Pathology), we demonstrate spatial imaging of neural cell and neurodegenerative disease-specific targets (beta-amyloid and phospho-tau species), and key post-translational modifications. Furthermore, advanced neuronal segmentation algorithms allow for specific tracing and segmentation of single neurons, including axons, across hundreds of micrometers of space. Ultimately, the ability to image and evaluate 5-fold larger sections of human brain tissues with unparalleled spatial context and high-plex analyte contents allows the grander scale of disease anatomy and processes to be cataloged with exquisite precision and accuracy.

Methods


Taking advantage of oligonucleotide barcode-conjugated antibodies and large-format flow cell, the NanoString® CosMx™ Spatial Molecular Imager (SMI) platform allows high-plex profiling of protein expression in FFPE human brain section in spatial context, capturing the single-cell variance over large spatial dimension.

CosMx 68-plex protein panel




Large format flow cell

- Imageable area > 1600 mm², >5x increase from the imageable area of standard CosMx SMI flow cell.
- Fluidic system on the CosMx SMI prototype was modified to enable uniform reagent exchange through the flow cell with minimal use of expensive reagents
- The system was able to achieve uniform reagent exchange with a reagent volume of only 250µL, despite the enormous imaging area of the flow cell




Sample preparation, flow cell assembly and target detection



High-plex Proteomics Profiles of Large Area Alzheimer's Disease Tissue



Alzheimer's disease (AD) is characterized by the spatiotemporal progressive deposition of pathologic peptides amyloid beta (Aβ) and hyper-phosphorylated tau (pTau). The burden of pTau, in the form of neurofibrillary degeneration, corresponds strongly with cognitive impairment and progresses generally from medial temporal lobe structures to neocortex. The middle temporal and superior temporal gyri (MTG and STG, respectively) represent a critical transition zone in disease progression, separating pTau accumulation associated with aging and pre-clinical AD from advanced stages of dementia-associated pTau pathology.

Spreading of pathology in AD




The novel Large Surface Area flow cell of the CosMx SMI prototype has > 1600 mm² of imaging area, which is a > 5-fold increase of standard CosMx SMI flow cell, and thus allows the study of pTau pathology within spatial context of both immediate and far-ranging neuroanatomy. Combining the large-format flow cell with a 68-plex CosMx SMI Human Neuroscience protein panel, we demonstrate spatial imaging of brain cell typing (GFAP, Iba1, NeuN), disease-specific targets (APP, Aβs, pTau), and key post translational modifications (pTau at multiple sites), tracking the spread of AD pathology across brain structures.

Brain section #1 from AD patient at Braak stage VI



Heterogeneity of plaque & tangle composition across brain regions

Brain section #2 from AD patient at Braak stage V



Conclusions

The large format CoxMx SMI technology provides the unique opportunity to study vulnerable cells continuously across multiple brain structures to understand the earliest neurodegenerative changes in disease. Using this approach, we have begun to shed light on unanswered questions such as how pTau species, comorbid protein-autopathies (such as TDP-43), and immune cells differentially impact vulnerable cell types and change dynamically change with disease progression.

By applying the 68-plex protein assay to AD tissue we observed:

- pTau pattern that tracks known spatiotemporal progression of tau
- Intriguing heterogeneity in plaque composition
- Differential subcellular location of pTau species


Critically, by sampling larger brain areas spanning multiple disease-associated regions from a single donor, we can begin to uncover the cellular and protein-level underlying key transition points along each stage of disease progression.

References


- He, S., Bhatt, R., Brown, C., et al. (2022). High-plex imaging of RNA and proteins at subcellular resolution in fixed tissue by spatial molecular imaging. *Nature Biotechnology*. 40(12):1794-1806. doi:10.1038/s41587-022-01483-z
- Stringer, C., Wang, T., Michaelos, M., & Pachitariu, M. (2021). Cellpose: a generalist algorithm for cellular segmentation. *Nature methods*, 18(1), 100-106.
- W. Tsai, "Moment-preserving thresholding: a new approach," *Computer Vision, Graphics, and Image Processing*, vol. 29, pp. 377-393, 1985.
- Braak H, Rüb U, Schultz C, Del Tredici K. Vulnerability of cortical neurons to Alzheimer's and Parkinson's diseases. *J Alzheimers Dis*. 2006;9(3 Suppl):35-44. doi:10.3233/jad-2006-9s305
- Brettschneider J, Del Tredici K, Lee VM, Trojanowski JQ. Spreading of pathology in neurodegenerative diseases: a focus on human studies. *Nat Rev Neurosci*. 2015;16(2):109-120. doi:10.1038/nrn3887
- Ahmadian SS, Rezvanian A, Peterson M, et al. Loss of calbindin-D28k is associated with the full range of tangle pathology within basal forebrain cholinergic neurons in Alzheimer's disease. *Neurobiol Aging*. 2015;36(12):3163-3170. doi:10.1016/j.neurobiolaging.2015.09.001

Robust Cell Segmentation Pipeline using High-plex multi-omics Images

(I) Preprocess Protein & RNA images for multi-omics assays.




(II) Perform high-plex segmentation leveraging protein colocalization & ML models.



Improved cell segmentation outcomes

- Green outline - Previous segmentation pipeline
- Red outline - New pipeline using high-plex protein images
- Blue shading - DAPI
- Gray shading - Max projection of high-plex protein images



Scan here to download or learn more



The CosMx™ SMI and decoder probes are not offered and/or delivered to the following UPC member states* for use in these countries for the detection of RNA in a method used for the detection of a plurality of analytes in a cell or tissue sample without the consent of the President and Fellows of Harvard College (Harvard Corporation) as owner of the Unitary Patent EP 4 108 782 B1. The use for the detection of RNA is prohibited without the consent of the of the President and Fellows of Harvard College (Harvard Corporation) as owner of the German part of EP 2 794 928 B1. The use for the detection of cellular RNA, messenger RNA, microRNA, ribosomal RNA and any combinations thereof in a method used in fluorescence in situ hybridization for detecting a plurality of analytes in a sample without the consent of the President and Fellows of Harvard College (Harvard Corporation) as owner of the German part of EP 2 794 928 B1. The use for the detection of cellular RNA, messenger RNA, microRNA, ribosomal RNA and any combinations thereof is prohibited without the consent of the of the President and Fellows of Harvard College (Harvard Corporation).