


Disentangling the relationship between pathology and cell-specific vulnerability using large format highly multiplexed in situ single-cell profiling of human brain

Victoria M. Rachleff¹, Emily E. Killingbeck Schneidereit¹, Chi Phan², Alex Bosworth², Arya Bahrami², Aster Wardhani², Lidan Wu², Mithra Korukonda², Ashley Heck², Claire Williams², Alyssa Rosenbloom², Joseph M. Beechem², and C. Dirk Keene¹


¹University of Washington Department of Laboratory Medicine and Pathology. ²NanoString® Technologies, a Bruker Company. Seattle, WA, USA.

BACKGROUND




Transition zones in AD

- Allocortical to neocortical progression
- Entorhinal cortex
- Hippocampus
- Temporal cortex (inferior → middle → superior temporal gyri)
- Frontal cortex



LARGE SURFACE AREA FLOW CELL


The CosMx™ SMI Spatial platform efficiently handles highly multiplexed protein analysis on FFPE tissues. Here we demonstrate a novel Large Surface Area flow cell with > 1600 mm² of imaging area, a > 5-fold increase. Combining the Large Surface Area flow cell imaging capabilities with a human specific > 68 plex CosMx™ Neural Cell Profiling and Alzheimer's Pathology Protein Panel, we demonstrate spatial imaging of brain cell typing (GFAP, Iba1, NeuN), disease-specific targets (APP, Aβ, pTau), and key post translational modifications (pTau at multiple sites), tracking the spread of AD pathology across brain structures.




- Fluidic system on the CosMx™ prototype was modified to enable uniform reagent exchange through the flow cell with minimal use of expensive reagents
- The system was able to achieve uniform reagent exchange with a reagent volume of only 250μL, despite the enormous imaging area of the flow cell

AD SPATIAL PROGRESSION

Several key brain regions were selected to capture the continuous progression of pathological peptides within a single brain donor. Regions 1 and 2 represent the regions earliest affected by tau pathology, while regions 3 and 4 represent cortical regions affected later in disease progression.




68-PLEX PROTEIN PANEL




Cell types and pathological peptides each have their own distinct roles in the pathophysiological progression of AD. Glia, such as astrocytes (GFAP) and microglia (IBA1), increase their abundance and aggregate in areas of high pathological peptide burden, whereas neurons (MAP2), particularly excitatory neurons, appear to decrease in relative abundance with increasing pathological stage. The NanoString CosMx high plex protein panel enables the assessment of cell subtypes as well as intracellular (such as pTau) and extracellular (such as Aβ) pathological peptides. Further, the panel captures a range of phosphorylation states for key AD proteins like pTau.


METHODS



Sample preparation, flow cell assembly and target detection



CELL TYPING AND SEGMENTATION



FUTURE DIRECTIONS

Three core components are essential to disentangle the relationship between pathology and cell-specific vulnerability: spatial context, pathological protein distribution, and highly resolved cell type labels. This dataset provides a detailed yet large-scale protein landscape. Combining this dataset with NanoString's 6K RNA panel will unlock the full potential of the technology and complete this goal.

FOVs (highlighted in red) selected for profiling on the NanoString 6K RNA panel that is currently being run on the same tissue section.

REFERENCES

[1] He, S., Bhatt, R., Brown, C., et al. (2022). High-plex imaging of RNA and proteins at subcellular resolution in fixed tissue by spatial molecular imaging. *Nature Biotechnology*. 40(12):1794-1806. doi:10.1038/s41587-022-01483-z

[2] Stringer, C., Wang, T., Michaelos, M., & Pachitariu, M. (2021). Cellpose: a generalist algorithm for cellular segmentation. *Nature Methods*, 18(1), 100-106.