



Highest plex spatial multiomics as a discovery tool for senescence cell biology in Alzheimer's disease




Hannah R. Hudson¹, Timothy C. Orr², Alyssa Rosenbloom³, C. Dirk Keene⁴, Joseph Beechem³, Miranda E. Orr^{2,5}

¹Translational Neuroscience, Wake Forest University School of Medicine, Winston-Salem, NC, ²Internal Medicine Geriatrics and Gerontology, Wake Forest University School of Medicine, Winston-Salem, NC, ³NanoString Technologies (a Bruker company), Seattle, WA, ⁴Lab. Medicine and Pathology, Univ. of Washington, Seattle, WA, ⁵Salisbury VA Medical Center, Salisbury, NC

Abstract

- Intraneuronal deposition of aberrant tau protein into neurofibrillary tangles (NFTs) closely correlates with neurodegeneration and cognitive decline. Pathogenic tau induces cell death in some neurons while other neurons evade apoptosis by entering cellular senescence, a complex stress response leading to a change in cell fate¹⁻⁴. We refer to senescent neurons as "neurocent" and are working to understand the role of tau in driving this cell fate.
- Removing senescent cells improves brain structure and function underscoring the importance of exploring senescent cells, their consequences across the brain and potential molecular targets for interventions¹⁻⁴.
- The advancement of spatially resolved, multiplex proteomic and transcriptomic technologies has revolutionized and redefined the approaches to complex biological questions pertaining to tissue heterogeneity, microenvironments, cellular interactions, cellular diversity, and therapeutic response in neurodegenerative diseases involving tau^{5,6}.
- Spatial transcriptomics has traditionally led the way in plex, multiple studies have demonstrated a poor correlation between RNA expression and protein abundance, owing to transcriptional and translational regulation, target turnover, and most critically, post-translational protein modifications. Therefore, a more holistic, ultra-high-plex, and high-throughput proteomic atlas approach becomes critical for the next phase of discovery biology.
- A Digital Spatial Profiler platform is uniquely suited to support high-plex proteomics, allowing for the simultaneous analyses of proteins from discrete regions of interest (ROIs) in FFPE tissue sections while preserving spatial context.
- Here, we use a barrier-breaking, highest plex spatial multiomics solution to explore the connection between tau neuropathology and senescence in human Alzheimer's disease. We also profiled human brains from patients with amyotrophic lateral sclerosis (ALS) and healthy controls to help identify disease-specific profiles.


Profiling Alzheimer's disease pathologies across the primary motor cortex and visual cortex




Orr Custom 63-Plex

| Phospho-Tdp-43 (S409+S410) | RIPK1 | GZMA | Phospho-RIPK3 (S164 + T165) |
|----------------------------|---------------------------|--------------------------|--------------------------------|
| DAT | cGAS | DE | Par7 |
| 4R tau | Phospho-Tau (S262+T263) | LRRK2 | PINK1 |
| Phospho-Tau (T181) | IBA1 | MAP2 | PLCG1 |
| Phospho-MLKL (S358) | Phospho-PRAS40 (T246) | Phospho-MEK1 (S217+S221) | PSEN1 |
| Map6 | ApoA1 | Phospho-Tau (S214) | NeuN |
| FIP200/191ccc1 | APOE | Phospho-Tau (S199) | Phospho-Alpha-synuclein (S129) |
| P2RX7 | Amyloid Beta 1-40 | Phospho-Tau (T231) | Alpha-synuclein |
| Cleaved Caspase 9 | Amyloid Precursor Protein | Phospho-Tau (S396) | STING |
| ParK2 | Amyloid Beta 1-42 | Phospho-Tau (S404) | Synaptophysin |
| RIPK3 | BACE1 | Tau | Tdp-43 |
| Phospho-RIPK1 (S161) | BAG3 | Myelin basic protein | TFEB |
| TREM2 | Calbindin | Nephrin | TMEM119 |
| SIRP-alpha | FUS | Neurofilament light | ParK5 |
| MLKL | GFAP | NRGN | VP35 |
| GBA | GFPMB | Olfg2 | |


Ultra-high-plex multiomics reveals distinct microenvironments immediately surrounding plaques specific to disease and brain region




Single cell ultra-high-plex proteomics demonstrates distinct tau tangle profiles within cortical layers in the primary motor cortex of human Alzheimer's disease





GeoMx Digital Spatial Profiler



Highest Plex Spatial Multiomics



Spatial Proteomic Workflow #2: Seq Code RNA Assay and Pro Code Protein Assay



Summary

- GeoMx DSP Platform Capabilities:**
- Single neuron resolution of multiple analytes across entire FFPE tissue sections.
 - Spatially-resolved multiomic data including whole transcriptome + 643-plex proteins with or without post-translational modifications.
- Preliminary Findings and On-Going Analyses:**
- Neuropathologies (NFTs and Aβ plaques) display molecular signatures that differ across diseases (AD with or without neuritic plaques, ALS and control), between cortical regions, and between cortical layers within the same region.
 - Whole transcriptome analyses reveals spatially resolved, co-localized gene modules in the plaque environment that differ across diseases and brain regions.
 - Evaluating RNA and protein expression in single neurons across pathologies suggests low concordance; analyses are ongoing.
 - Combining IPA with the neuro panel is providing insights into neuroscience at an unprecedented molecular detail.

References

- Musi N et al. *Aging Cell*. 2018.
- Dehtori SK et al. *Nat Aging*. 2021.
- Gonzales MM et al. *Nat Med*. 2023.
- Riessland M et al. *Nat Rev Drug Discov*. 2024.
- Walker J et al. *Acta Neuropathol Commun*. 2022.
- Walker J et al. *Alzheimers Dement*. 2024.

Funding (to MEO): Department of Veterans Affairs I01BX005717, NIA R01AG068293, Cure Alzheimer's Fund, AFAR/Hevolution, Rainwater Charitable Foundation).