



# Highly multiplexed spatial single-cell multiomic imaging of tau neuropathology in human brain tissue sections




Timothy Orr<sup>1</sup>, Alyssa Rosenbloom<sup>2</sup>, Ashley (THE DRAGON) Heck<sup>2</sup>, Lidan Wu<sup>2</sup>, C.S. Latimer<sup>3</sup>, C. Dirk Keene<sup>3</sup>, Joseph Beechem<sup>2</sup>, \*Miranda Orr<sup>1,4</sup>

<sup>1</sup>Wake Forest Univ. Sch. of Med., Winston-Salem, NC; <sup>2</sup>NanoString Technologies, a Bruker company, Seattle, WA  
<sup>3</sup>Lab. Medicine and Pathology, Univ. of Washington, Seattle, WA; <sup>4</sup>Salisbury VA Medical Center, Salisbury, NC


## Introduction

- Intraneuronal deposition of aberrant tau protein into neurofibrillary tangles (NFTs) closely correlates with neurodegeneration and cognitive decline. Pathogenic tau induces cell death in some neurons while other neurons evade apoptosis by entering cellular senescence, a complex stress response leading to a change in cell fate<sup>1,2</sup>. Removing senescent cells improves brain structure and function underscoring the importance of exploring senescent cells, their consequences across the brain and potential molecular targets for interventions<sup>1,3,4</sup>.
- To study the transition of post-mitotic neurons into senescence, along with its impact on neighboring cells and tissue-wide effects, high-resolution spatial analyses are essential<sup>5,6</sup>. These analyses should capture changes in protein levels, post-translational modifications (e.g., phosphorylation), and gene expression. A spatial molecular imager (SMI) provides ultra-high plex detection of RNA and protein, enabling detailed molecular profiling within a spatial context.
- Here we showcase the multi-omic capabilities of the SMI in detecting the entire transcriptome and >68 proteins from healthy and diseased human brain sections.




## High-plex spatially resolved protein data enables analyses of pathogenic tau and its effects on neighboring cells


### Unsupervised clustering: 35 clusters encompassing major neural cell types




### Clustering on neighborhood cell type composition identified tissue niches matching anatomical structures and enabling DE analyses across regions




### Differential pTau between diseases




### Spatial patterns revealed for distinct pTau species in ALS




### Protein object analyses: AT8 segmentation and differential expression at distances from NFTs



## SMI protein and RNA data reveal G<sub>x</sub> neuroscence




Previously discovered neurescent genes spatially map and colocalize with NFTs. Left: Previous studies developed senescence eigengenes to identify senescent cells in single nucleus RNA sequencing datasets. The eigengenes indicated that NFTs-bearing neurons were a primary source of senescent cells in the AD brain. The genes within the eigengene shown to contribute most weight to the senescence phenotype included *ETS2*, *RBL2*, *E2F3*, *ATM* and *CDK6*. CosMx SMI indicate that these genes are upregulated in NFT-bearing neurons. Further differential expression analyses identified a novel gene, *UCCRHL*, that also spatially maps and co-localizes with NFTs (green). Right: Differential gene expression analyses of CosMx SMI datasets has revealed additional, novel, putative neuronal senescence, "neurescence" genes.





Cell cycle events and cell fate options of postmitotic neurons. Left: Neurons differentiate from precursor cells, exit the cell cycle and may 1) remain in G<sub>0</sub> through effective cell maintenance, damage repair and suppressing cell cycle events; 2) undergo cell death; 3) re-engage the cell cycle where they may undergo cell death, continue through various stages of the cell cycle including G<sub>1</sub> and S or; 4) enter senescence (G<sub>0</sub>). Entry to G<sub>0</sub> may occur directly from G<sub>0</sub> (4a) or by exiting the cell cycle after stress-induced re-entry (4b-4d). CDK1: cyclin dependent kinase inhibitor; SASP: senescence associated secretory phenotype; NASP: neurescence associated secretory phenotype; \*Cell death: multiple neuronal cell death mechanisms, beyond apoptosis, have been described for neurodegenerative diseases. Right: Neurons with NFTs (red) co-express the G<sub>0</sub> gene, *CDKN2D* (white dots), which encodes p19<sup>INK4d</sup>. Green: GFAP (astrocytes); Blue: DNA (nuclei); Gray: Histone (nuclei).

## CosMx™ SMI Protein & RNA Workflows




### Differential protein expression in entorhinal cortex neurons




## Whole transcriptome RNA data identifies spatially colocalized gene modules that differentiate disease from control tissue

### Unsupervised clustering identified 22 clusters, encompassing major neural cell types



### InSituCor identifies spatially colocalized gene modules that differ with disease



UMAP visualization showing 22 distinct clusters representing major neural cell types from tissues affected by amyotrophic lateral sclerosis (ALS), Alzheimer's disease (AD) and control tissues. Control tissues are from two different ages. Young: 36 years-old and Old: 73 years-old. Clusters reflect similarities and differences in cell type profiles across disease states and age groups. Spatial plots depict the tissue distribution of the cell types identified in the UMAP analysis, showing their physical locations within the tissue sections.

InSituCor detects gene expression patterns that are spatially related, while filtering out trivial correlations, such as those that might occur simply due to cell type similarities. Applying InSituCor to the CosMx SMI dataset revealed multiple gene modules that were unique to (a) ALS entorhinal cortex; (b) young control brain; (c-d) AD brain with either elevated as shown in (c) or decreased expression as shown in (d); (e) control; and (f) ALS brain.

## Summary

- CosMx SMI Platform Capabilities:**
- 50nm resolution of analytes across entire FFPE tissue sections
  - Spatially-resolved multiomic data including whole transcriptome + 68-plex protein, and post-translational modifications
  - Identification of major brain cell types and spatial niches
  - Whole transcriptome analyses reveals spatially resolved, co-localized gene modules
  - Protein object analyses reveals molecular, cellular and niche changes in relation to neuropathology
- Preliminary Findings and On-Going Analyses:**
- Molecular signatures of distinct neuropathologies (NFTs and Aβ plaques)
  - Spatial resolution and intersection of neuropathology and neurescent cells
  - Multi-analyte, cell and niche specific changes and vulnerabilities to distinct diseases (ALS vs AD)

## References

- Musi N et al., *Aging Cell*. 2018.
- Dehkordi SK et al., *Nat Aging*. 2021.
- Gonzales MM et al., *Nat Med*. 2023.
- Riessland M et al., *Nat Rev Drug Discov*. 2024.
- Walker J et al., *Acta Neuropathol Commun*. 2022.
- Walker J et al., *Alzheimers Dement*. 2024.