


# In situ 3D spatial characterization of genome architecture and whole transcriptome imaging in single cells with the PaintScape™ system and CosMx® SMI on ER+ breast cancer cells



Huy Nguyen<sup>1</sup>, Yi Cui<sup>2</sup>, Brian Smart<sup>1</sup>, Austin Nguyen<sup>1</sup>, Andrea Floris<sup>1</sup>, Serdar Tulu<sup>1</sup>, David Castillo<sup>1</sup>, Sierra McKinzie<sup>2</sup>, Isabel Lee<sup>2</sup>, Shanshan He<sup>2</sup>, Shyamtanu Chatteraj<sup>1</sup>, Jude Dunne<sup>1</sup> and Joseph Beechem<sup>2</sup>  
<sup>1</sup>Bruker Spatial Genomics Inc., USA; <sup>2</sup>Bruker Spatial Biology, USA



- Here we present a novel jebFISH™ protocol on the PaintScape™ platform and the CosMx® Spatial Molecular Imager (SMI) Whole Transcriptome (WTX) assay used to characterize differences in genome structure, gene expression and phenotype at subcellular resolution between ER+ breast cancer cell line MCF7 and ER- normal breast cell line MCF10a.
- ER+ cell lines often exhibit selective amplification of chromosomal regions (e.g. 17q23 and 20q23 in MCF7) containing distant estrogen response elements (DEREs), which facilitate long range chromatin interactions causing transcriptional repression of tumor-suppressor genes and activation of oncogenes.<sup>1</sup>
- ER activation induces upregulation in key proliferative and differentiation gene targets, including CCND1, MYC through both proximal ERs and long-range DEREs that coordinate enhancer-promoter interaction.<sup>2-3</sup>
- Luminal epithelial and non-invasive identity of MCF7 is maintained due to high expression of GATA3, a core luminal transcription factor co-expressed and co-regulated with ER.<sup>4</sup>
- The ER-GATA3 axis maintains luminal epithelial programs by maintaining high CDH1 expression while transcriptionally suppressing EMT drivers such as SNAI1 and SNAI2, despite higher genomic copy number of these two genes in MCF7.<sup>5-6</sup>
- Loss of ER/GATA3 signaling is associated with EMT activation and more aggressive, less differentiated tumor phenotypes<sup>7</sup>
- Selective copy gain, rearrangement, and altered 3D architecture of those oncogenes and their regulatory elements cause dysregulation of oncogenic signaling pathways and increases cancer cell proliferation providing positive selection for cancer progression.

## Technology and Methods

**A The PaintScape System powered with jebFISH™ technology to visualize in-situ 3D genome organization in situ in single cells**

**What is the PaintScape System?**

- An instrument platform that combines state-of-the-art optics, sophisticated bioinformatics, and a reagent system based on powerful **jebFISH** chemistry for direct *in situ* visualization of the 3D genome in individual cells

**What is jebFISH?**

- High-plex, efficient** chemistry for *in situ* visualization of the 3D genome in single cells (in cell lines and fresh frozen tissue) using a proprietary multiplex optical **barcoded** chemistry
- Proprietary **jebSmart™** (smart barcode) method to allow precision loci identification and localization

**B OncoPaint™ HuCL Cancer Panel Kit**

Each gene region is painted with 15 loci. The MDM2 gene region is shown as an example.

**>1000 plex jebFISH painting of MCF10a**

PaintScape tracing of 1050 target loci in MCF10a. Each chromosome is painted as described in methods, section B


**OncoPaint™ HuCL Panels Used: Transcriptional Regulation, DNA Repair, Cell Cycle and Apoptosis, Chromatin Structure and A/B ploidy panels**

**C CosMx WTX assay design**



CosMx RNA detection schematic

WTX chromosome coverage


WTX RNA detection in cultured MCF7 cells




## PaintScape system characterizes and directly visualizes 3D genome organization and structural genotype of ER+ MCF7 cells



## PaintScape 3D genome structural data yields insights into differential gene expression as measured by CosMx WTX



## PaintScape 3D genomics and CosMx WTX characterizes ecDNA element driven oncogenic expression in MCF7



## Conclusions

- The PaintScape system enables high-resolution identification, characterization, and direct *in situ* visualization of single-cell 3D genome architecture in single cells. CosMx WTX profiles spatial gene expression at a single-cell level revealing disease state specific oncogenic function. Here, the combined approach of PaintScape and CosMx reveals how disruptions in genome organization drive disease mechanisms in ER+ breast cancer by linking altered chromatin topology to dysregulated expression of key oncogenic pathway genes. Distinct structure-function correlations are revealed:
  - Copy number variation including amplification and deletions *in situ* in single cells and sub-population
    - Specific gain in Chr 17q and 20q regions containing DEREs in MCF7 BC cells
    - Focal amplifications of important oncogenes BRIP1, MYC, SNAI1
    - Loss in CDKN2A gene on Chr 9
  - Single cell ploidy variation of whole chromosomes and sub-chromosomal regions
    - Higher variation in sub-population of cells with variable chromosome ploidy in cancer cells compared to normal cells
    - Chr17 and Chr20 show much higher chromosome ploidy in cancer cells
  - Inter-chromosomal interactions, potential translocations and simultaneous proximity of multi-loci interactions *in situ* in single cells
    - Close proximity of chromosomes 3, 5 and 9 regions in MCF10a consistent with known translocations t(3;5) and t(3;9)
    - Close proximity of GATA3 (Chr 10) and HOXA (Chr 7) loci in MCF7 causing dysregulated gene expression of GATA3 and certain HOXA genes in MCF7
- ecDNA copy number, structure, nuclear locations and *in situ* interactions
  - A fraction of amplified DERE targets of Chr17q and Chr20q regions were further from Chr17 or Chr20 territory suggesting possible ecDNA nature of those amplified targets
  - BRIP1 gene region shows very high copy number due to ecDNA amplification in MCF7
  - A significant fraction of the BRIP1 ecDNA formed self-interacting hub like structure in single MCF7 cells
  - BRIP1 ecDNA closely interacts with amplified DEREs on Chr20q including SNAI1 oncogene
- 3D genome structure, spatial gene expression and multiomic structure-function correlations
  - CosMx WTX identified elevated expression of ER-driven oncogenes GATA3 and CCND1 in MCF7 cells, consistent with enhanced estrogen-dependent luminal signaling and proliferative activity
  - PaintScape 3D genome data shows increased enhancer-promoter interactions and strengthen intra-TAD interactions for these genes revealing mechanism of higher expression
  - PaintScape system characterized and directly visualized unique ecDNA structure containing relevant oncogenes BRIP1, TBX4, MED13 in MCF7
  - CosMx WTX shows higher expression of these genes in MCF7 compared to MCF10a revealing functional relevance of ecDNA amplification which promotes opportunistic enhancer hijacking within ecDNA elements and promoting dysregulated expression

## PaintScape™ System and CosMx WTX enables *in situ* direct visualization of 3D genome changes and spatial gene expression in single breast cancer cells

### References

- Cancer Res 2009;69(14):5946-53
- Mol Cell Biol 2004, 24, 7260-7274
- Mol Endocrinol 2011, 25, 1527-1538
- Cancer Res 2007, 67, 6477-6483
- Oncogene 2010, 29, 1451-1462
- Cancer Advances 2025, 8, e25016
- J Biol Chem 2010, 285, 14042-14051
- Nat Commun 2024, 15, 6130
- GeneCards.org
- Nature 2021, 600, 731



Scan to download or learn more