



Integrated 3D genome, transcriptome and proteomic assays – simultaneously observing genome architecture and transcriptional activity at single-cell resolution





Jude Dunne¹, David King¹, Nagasangeethavani Anandakrishnan¹, Johnson Huynh¹, Serdar Tulu¹, Wyatt Eng¹, Marc Glazer¹, Brian Smart¹, David Castillo¹, Kenny Chung¹, Huy Nguyen¹, Shyamtanu Chatteraj¹, and Doug Werner¹

¹Braker Spatial Genomics Inc., Research and Development, San Jose, California, USA





- Here we present a novel jebFISH™ protocol on the PaintScope platform used to characterize differences in genome structure, gene expression, and phenotype at subcellular resolution between HER2⁺ breast cancer cell lines HCC1954 and SKBR3.
- Variation in chromatin structure, including whole or partial chromosome copy gain or loss, genome openness or compartment switching, and TAD boundary disruptions can directly alter transcriptional activity and protein expression.
- Cancer cells heavily rely on dysregulated cell cycle and signaling pathways for proliferation and invasion to further fuel abnormal cell growth and epithelial-mesenchymal transition. These signaling pathways consist of a complex network of interactions involving several key oncogenes.
- HCC1954 and SKBR3, both HER2⁺ breast cancer cell lines used in this study, contain hotspots of clustered breakpoints on specific chromosomal DNA or extra chromosomal DNA (ecDNA) regions and often show complex focal amplifications¹⁻⁴. Selective copy gain, rearrangement, and altered 3D architecture of underlying oncogenes and their regulatory elements cause dysregulation of oncogenic signaling pathways and increases cancer cell proliferation providing positive selection for cancer progression.
- We present a multiomics assay on the PaintScope System using our OncoPaint™ Oncogenic Pathways Panels with associated transcriptome targets that seeks to characterize these structural genomic changes and the impact they have on gene expression.
- Automated imaging decodes both DNA and RNA targets *in situ*. Optical barcodes label more than 1,000 genomic loci across the genome, including dense coverage of 45 gene regions spanning multiple cancer-relevant signaling pathways. The panel also includes probes for transcripts from key cancer-associated pathways, enabling simultaneous assessment of pathway-specific transcriptional outputs. Together, this integrated approach allows analysis of chromatin structure and gene expression within the same cells, enabling direct comparison between local chromatin organization and transcriptional activity.



Technology and Methods




PaintScope system directly visualizes 3D genome structure copy number and interactions between genes *in situ* in single breast cancer cells



PaintScope platform correlates transcriptomics with 3D genome structure in different breast cancer cells



PaintScope system enables assessment of genome structure in cell subpopulation based on transcriptomic expression analysis of same cells



Conclusions

- The PaintScope system enables high-resolution identification, characterization, and direct *in situ* visualization of single-cell 3D genome architecture and spatial gene expression in single cells revealing disease state specific oncogenic function. The novel same-cell multiomic approach on the PaintScope platform shows how disruptions in genome organization drive disease mechanisms in HER2⁺ breast cancer by linking altered chromatin topology to dysregulated expression of key oncogenic pathway genes. Distinct structure-function correlations are revealed.
- 3D genome structure, spatial gene expression and same cell multiomic structure-function correlations *in situ*
- PaintScope spatial transcriptomics identified differential expression of key breast cancer related genes CDH1, SOX9, STAT3, FOXM1 that are associated with cell state, stemness and oncogenic transcriptional reprogramming between HCC1954 and SKBR3
- This is driven by underlying copy number variations (CNVs), genomic rearrangements, structural alterations
- Copy number variation including amplification and deletions *in situ* in single cells and sub-population
 - Focal amplifications of important oncogenes: MYC, CCND1, SNAI1 in HCC1954; MET, MYC and SNAI1 in SKBR3
 - Loss in CDH1 DNA copy in SKBR3
- Single cell ploidy variation of whole chromosomes and sub-chromosomal regions
 - Higher variation in sub-population of cells with variable chromosome ploidy in HCC1954 than SKBR3
- ecDNA copy number, structure, nuclear locations and *in situ* interactions
 - ecDNA type amplification of MYC and CCND1 genes in HCC1954 cells; MET, MYC and SNAI1 genes in SKBR3 cells
 - These ecDNAs show simultaneous proximity of multi-loci interactions involving key oncogenes in individual single B0 cell sub-populations
- Same cell multiomic ability of PaintScope system identified sub-population of cancer cells with distinct 3D genome structure and spatial RNA expression profile
 - PaintScope transcriptomics identified high and low FOXM1 RNA expressing cells in SKBR3
 - High FOXM1 expressing cells show well defined TAD structure with increased enhancer-promoter interactions
 - These findings establish a framework for linking spatial transcriptomic signatures to structural genomic events, with implications for understanding intratumoral heterogeneity and clonal evolution in breast cancer cell lines

PaintScope System directly visualizes *in situ* 3D genome changes and associated RNA expression in the same cells

References

- Cancer Genetics 204 (2011) 447-457
- Genome Res. 28 (2018) 1126-1135
- Genome Biol 20 (2019) 27
- bioRxiv preprint doi: https://doi.org/10.64898/2025.12.01.691649
- Nature 635 (2024) 193-200
- Nat Genet 55 (2023) 2189-2199
- Nature 600 (2021) 731-736
- Genes Chrom Cancer 46 (2007) 543-558
- Life Sci Alliance. 21 (2023) e202302154
- Int. J. Oncol 35 (2009) 57-68