


# PaintScape™ enables multiomic in-situ direct visualization of spatial 3D genome architecture of single cells in intact fresh frozen Glioblastoma tissue in native tumor microenvironments


Huy Nguyen<sup>1</sup>, Matija Snuderl<sup>2</sup>, Yuxiu Wang<sup>2</sup>, Neha Dhasmana<sup>2</sup>, Pamela Flatley<sup>1</sup>, Stacy Elliott<sup>1</sup>, Doug Werner<sup>1</sup>, Jude Dunne<sup>1</sup> and Shyamtanu Chatteraj<sup>1</sup>  
<sup>1</sup>Brucker Spatial Genomics Inc., USA, <sup>2</sup>Department of Pathology, NYU Grossman School of Medicine, New York, New York, USA

## Introduction




- Here we present a novel jebFISH™ tissue genome painting protocol on the PaintScape™ platform to reveal the 3D genome organization *in situ* of fresh frozen (FF) glioblastoma (GBM) tissue at single cell, sub-population, and population levels.
- GBM is the most aggressive and common type of brain tumor and is characterized by high intra- and inter-tumoral molecular heterogeneity. Over 70% of GBM patients exhibit genetic signatures, such as Chr7 amplification and Chr10 loss of heterozygosity (LOH)<sup>1</sup>. Also, GBM patient outcome is strongly impacted by tumor microenvironments, with hypoxia noted as a significant detrimental factor<sup>2</sup>.
- Chromosome conformation capture analysis on bulk tissue has identified the role of the 3D genome in GBM<sup>3</sup>, but further understanding of GBM requires tools that can study the genomic complexity with single cell resolution and spatial context.
- Here we use the PaintScape System to directly visualize the genome organization of GBM tissue *in situ* in single cells of a 30-µm FF tissue section. In addition, hypoxic and normoxic regions of the tissue section were identified to perform further structural genome organization analysis and comparison of tumor cells located in those regions and reveal differences or similarities in the respective structural genome landscapes of these sub-regions.
- The powerful nature of PaintScape-derived data enables even further filtering to focus on subpopulations of cells with specific genome organization traits for even deeper analysis of combinations of parameters and combinations of specific cell subsets that were not possible to compare before. For example, here, the 3D genome landscapes of a subset of tumor cells with Chr 7 trisomy from both hypoxic and normoxic regions of the GBM tissue is selected and further compared—revealing specific structural genome differences that contribute to the molecular heterogeneity of GBM tumor cells.


## Technology and Methods




## PaintScape system directly visualizes and characterizes *in situ*, single cell 3D genome organization and structural genotype of Glioblastoma from 30 µm thick fresh frozen tissue section




## PaintScape characterizes 3D genome and conformational changes in hypoxic and normoxic GBM tissue environments



## PaintScape characterizes MET gene region amplification and conformational changes in hypoxic and normoxic regions



## Single cell 3D genome conformational differences in Chr7 trisomy subpopulations in hypoxic and normoxic GBM tissue



## Conclusions

Distinct structural genotypes and genome conformations revealed by jebFISH on the PaintScape System:

- Canonical GBM genomic features including Chr10 loss and subpopulations
- Hypoxic vs. normoxic:
  - CNVs in GBM and hypoxia related genes and 3D genome differences (Chr10 further from other chromosomes in hypoxic)
  - Chr7 3N cells: CNVs in GBM and hypoxia related genes and 3D genome differences (Chr7 and Chr10 differential interactions)
- High levels of MET amplification forming ecDNA with neighboring genes CAPZA2 and ST7 in hypoxic and normoxic regions

PaintScape™ enables *in situ* direct visualization of 3D genome in FF Glioblastoma tissue in spatially resolved tumor microenvironments

## References

- Cancer Res. 2024 Oct 15;84(20):3464–3477
- Cancers (Basel). 2022 Feb 24;14(5):1176
- Sci. Adv. 2011, eadn2830(2025).
- Int J Oncol. 2022 Jan 10;60(2):15
- Int J of Rad Oncol, Volume 90, Issue 1, S39
- Cancer Cell, Volume 26, Issue 2, 288–300
- JAMA. 2009 Jul 15;302(3):276–89
- Front Mol Neurosci. 2017 Oct 17;10:324
- Front Mol Neurosci. 10:324.
- Genes Dev. 2000 Feb 15;14(4):391–396
- Cancer Discov. 2025 Jun 15;15(6):1141–1158



Scan to download or learn more