

Integrated spatial multiomic profiling of 3D genome architecture and transcriptome-wide gene expression in breast cancer models

Yi Cui², Huy Nguyen¹, Brian Smart¹, Austin Nguyen¹, Andrea Floris¹, Serdar Tulu¹, David Castillo¹, Sierra McKinzie², Isabel Lee², Shanshan He², Shyamtanu Chatteraj¹, Jude Dunne¹ and Joseph Beechem²
¹Bruker Spatial Genomics Inc., USA; ²Bruker Spatial Biology, USA

Introduction

Visualize and analyze the 3D genome *in situ* at single-cell resolution with jebFISH™ technology

Aggregate to form sub and whole population analyses

- Here we present a novel jebFISH™ protocol on the PaintScope™ platform and the CosMx® Spatial Molecular Imager (SMI) Whole Transcriptome (WTX) assay used to characterize differences in genome structure, gene expression and phenotype at subcellular resolution between ER+ breast cancer cell line MCF7 and ER- normal breast cell line MCF10a.
- ER+ cell lines often exhibit selective amplification of chromosomal regions (e.g. 17q23 and 20q23 in MCF7) containing distant estrogen response elements (DEREs), which facilitate long range chromatin interactions causing transcriptional repression of tumor-suppressor genes and activation of oncogenes.¹
- ER activation induces upregulation in key proliferative and differentiation gene targets, including CCND1, MYC through both proximal EREs and long-range DEREs that coordinate enhancer-promoter interaction.^{2,3}
- Luminal epithelial and non-invasive identity of MCF7 is maintained due to high expression of GATA3, a core luminal transcription factor co-expressed and co-regulated with ER.⁴
- The ER-GATA3 axis maintains luminal epithelial programs by maintaining high CDH1 expression while transcriptionally suppressing EMT drivers such as SNAI1 and SNAI2, despite higher genomic copy number of these two genes in MCF7.^{5,6}
- Loss of ER/GATA3 signaling is associated with EMT activation and more aggressive, less differentiated tumor phenotypes⁷
- Selective copy gain, rearrangement, and altered 3D architecture of those oncogenes and their regulatory elements cause dysregulation of oncogenic signaling pathways and increases cancer cell proliferation providing positive selection for cancer progression.

Technology and Methods

A The PaintScope System powered with jebFISH™ technology to visualize in-situ 3D genome organization in situ in single cells

What is the PaintScope System?

- An instrument platform that combines state-of-the-art optics, sophisticated bioinformatics, and a reagent system based on powerful jebFISH chemistry for direct *in situ* visualization of the 3D genome in individual cells

What is jebFISH?

- High-plex, efficient chemistry for *in situ* visualization of the 3D genome in single cells (in cell lines and fresh frozen tissue) using a proprietary multiplex optical barcoded chemistry
- Proprietary jebSmart™ (smart barcode) method to allow precision loci identification and localization

B OncoPaint™ HuCL Cancer Pathways Kit 1 Panels

Each gene region is painted with 15 loci. The MDM2 gene region shown as an example.

OncoPaint™ Panel Modules Used: Transcriptional Regulation, DNA Repair, Cell Cycle and Apoptosis, Chromatin Structure and A/B ploidy panels

C CosMx WTX assay design

CosMx RNA detection schematic, WTX chromosome coverage, WTX RNA detection in cultured MCF7 cells

PaintScope system characterizes and directly visualizes 3D genome organization and structural genotype of ER+ MCF7 cells

A PaintScope system identifies differential 3D genome amplification and conformation of oncogenic regions between normal and ER+ breast cancer cells

Distributions of chromosome counts per cell show the chromosomal copy gain in MCF7

Population-level distance heatmap shows close proximity of chromosomes 3, 5 and 9 regions consistent with known translocations t(3;5;9) and t(3;9)

Population-level distance heatmap shows close proximity of HOXA9 and GATA3 loci as well as Chr20q interactions with Chr1,8,11 and 17

B PaintScope system identifies cells with target amplification in context of chromosomal copy number

Select genes (GATA3, SNAI1 and BRIP1) showing target counts similar to, slightly in excess of and largely in excess of chromosome count.

C PaintScope system identifies in situ 3D genome conformation in ER+ breast cancer cells

PaintScope directly visualizes cells with a range of target amplifications, GATA3, SNAI1 and BRIP1. Note that GATA3 shows interaction with HOXA9 locus. Note that BRIP1 shows ecDNA like amplification relative to Chr17⁸

CosMx WTX assay characterizes unique gene expression profile of ER+ BC cell MCF7 in situ and identified key signaling pathways and unique cell states compared to normal breast cell MCF10a

A MCF7 and MCF10a cells show transcriptome-wide differences, with key Estrogen-Responsive Genes (ERGs) differentially expressed

B MCF7 cells possess ER-DRIVEN molecular markers with specific oncogenic functions

C CosMx WTX reveals inter- and intrapopulation expression heterogeneity

D ER+ MCF7 cells display spatial nuclear accumulation of ERG transcripts consistent with ligand-induced chromatin activation and transcriptional engagement

PaintScope 3D genome structural data yields insights into differential gene expression as measured by CosMx WTX

A PaintScope in situ visualization contextualizes differences in gene expression between cell lines at cell population level

Differentially expressed genes in MCF7 such as GATA3 and CCND1 can be understood in terms of copy number gain and promoter-enhancer interaction proximity.

Gene Locus Distance Maps showing Topologically associated domains (TADs)

SNAI1 and BRIP1 are amplified in MCF7 and do not show commensurate increase in gene expression.

Close interaction between Chr17-20q elements

PaintScope visualizes Chr 17q-20q interactions between DERE elements known to exhibit promoter-enhancer hijacking in ecDNA elements and known to dysregulate expression^{1,8}

PaintScope 3D genomics and CosMx WTX characterizes ecDNA element driven oncogenic expression in MCF7

A PaintScope system paints the amplified Chr17q region with higher genomic resolution around BRIP1

OncoPaint DNA Repair Pathway panel paints BRIP1 region with high fidelity and maps to portions of Chr17 in an inferred ecDNA structure for BRIP1⁸. Genes associated with ecDNA structure are assessed for copy number and expression

MCF7_BRIP1 loci on ecDNA structure

Single cell in situ visualization of ecDNA around BRIP1 locus on chr17 shows targets aggregated in a hub like manner⁹. Genes associated with this ecDNA such as TBX4 and INST2 do not show elevated expression despite high copy numbers indicating possible loss or hijacking of enhancer regions enabling dysregulated expression. Genes such as MED13 show elevated expression and may benefit from presence in ecDNA structures.

Conclusions

The PaintScope system enables high-resolution identification, characterization, and direct *in situ* visualization of single-cell 3D genome architecture in single cells. CosMx WTX profiles spatial gene expression at a single-cell level revealing disease state specific oncogenic function. Here, the combined approach of PaintScope and CosMx reveals how disruptions in genome organization drive disease mechanisms in ER+ breast cancer by linking altered chromatin topology to dysregulated expression of key oncogenic pathway genes. Distinct structure-function correlations are revealed:

- ecDNA copy number, structure, nuclear locations and *in situ* interactions
 - A fraction of amplified DERE targets of Chr17q and Chr20q regions were further from Chr17 or Chr20q territory suggesting possible ecDNA nature of those amplified targets
 - BRIP1 gene region shows very high copy number due to ecDNA amplification in MCF7
 - A significant fraction of the BRIP1 ecDNA formed self-interacting hub like structure in single MCF7 cells
 - BRIP1 ecDNA closely interacts with amplified DEREs on Chr20q including SNAI1 oncogene
- 3D genome structure, spatial gene expression and multiomic structure-function correlations
 - CosMx WTX identified elevated expression of ER-driven oncogenes GATA3 and CCND1 in MCF7 cells, consistent with enhanced estrogen-dependent luminal signaling and proliferative activity
 - PaintScope 3D genome data shows increased enhancer-promoter interactions and strengthened intra-TAD interactions for these genes revealing mechanism of higher expression
- PaintScope system characterized and directly visualized unique ecDNA structure containing relevant oncogenes BRIP1, TBX2, MED13 in MCF7
 - CosMx WTX shows higher expression of these genes in MCF7 compared to MCF10a revealing functional relevance of ecDNA amplification which promotes opportunistic enhancer hijacking within ecDNA elements and promoting dysregulated expression
- Inter-chromosomal interactions, potential translocations and simultaneous proximity of multi-loci interactions *in situ* in single cells
 - Close proximity of chromosomes 3, 5 and 9 regions in MCF10a consistent with known translocations t(3;5;9) and t(3;9)
 - Close proximity of GATA3 (Chr10) and HOXA (Chr7) loci in MCF7 causing dysregulated gene expression of GATA3 and certain HOXA genes in MCF7

PaintScope™ System and CosMx WTX enables in situ direct visualization of 3D genome changes and spatial gene expression in single breast cancer cells

References

- Cancer Res 2009;69(14):5946-53
- Mol Cell Biol 2004; 24: 7260-7274
- Mol Endocrinol 2011; 25, 1527-1538
- Cancer Res 2007; 67, 6477-6483
- Oncogene 2010; 29, 1451-1462
- Cancer Advances 2025; 8, e25016
- J Biol Chem 2010; 285, 14042-14051
- Nat Commun 2024; 15, 6130
- Genecards.org
- Nature 2021; 600, 731

Scan to download or learn more