

Cross-platform same slide multi-omics reveals Treg heterogeneity and links to spatial niche variability in pre-cancerous colonic inflammation

Daniel Jimenez-Sanchez¹, Matthew H. Ingalls¹, Sanghamithra Korukonda², Brian Lane¹, Isabella Peshek³, Patrick Danaher², Prajan Divakar², Oliver Braubach¹, Parambir S. Dulai³

¹Bruker Spatial Biology, St. Louis, MO, USA. ²Bruker Spatial Biology, Seattle, WA, USA. ³Division of Gastroenterology, Northwestern University, Chicago, Illinois, USA.

Introduction


FOXP3⁺ regulatory T cells (Tregs) maintain immune tolerance, and their dysfunction contributes to autoimmunity, chronic inflammation, and cancer. In ulcerative colitis, Treg heterogeneity is well documented and linked to colon cancer. However, Treg states are typically defined by transcriptional profiles, which may be confounded by discordance between RNA and protein levels.

Here, we use same-slide multi-omics with CellScape™ precise spatial proteomics and CosMx® Spatial Molecular Imager to quantify FOXP3 protein alongside gene expression. We find that proteomics reveals distinct FOXP3 protein levels defining Treg states within different spatial neighborhoods and functional contexts, which are further confirmed by spatial transcriptomics.


Methods

Eight FFPE intestinal biopsies from a pre-cancerous condition (ulcerative colitis) were profiled with CellScape™ (34-plex VistaPlex), followed by CosMx® Human Universal Cell Characterization Panel (1K-plex) on the same slide. Whole-slide multimodal alignment used shared morphology markers and WsiReg to generate affine registrations. Segmentation masks were harmonized for pixel-level correspondence, enabling direct integration of protein and RNA per cell. CD4⁺ T cells were classified into FOXP3_{low}, FOXP3_{mid}, FOXP3_{high} groups via quantitative automatic cell typing. Spatial neighborhoods (k=15) yielded 25 FOXP3-defined niches, and both protein and RNA features were included in neighborhood and state enrichment analyses.


FOXP3_{low} niches are enriched with stromal-vascular endothelial architecture




FOXP3_{mid} niches are enriched CD56+GranzymeB+ cells and remodeling vasculature



FOXP3_{high} niches are enriched in adaptive immune cells and proliferating T cells



FOXP3 expression levels map onto distinct spatial and phenotypic cellular neighborhoods




Same slide multi-omics captures protein and RNA from a single tissue section



Fig. 1 Same-slide CellScape (34-plex VistaPlex, protein) and CosMx (1K-plex, RNA) profiling on FFPE tissue enables direct, high-fidelity, molecular integration, eliminating the need for serial sections.

Shared nuclear and morphology markers enable pixel-level CellScape-CosMx registration



Pixel-level alignment enables shared:

- Cell Segmentation
- Cell Typing
- Protein Expression

Fig. 2 Whole-slide affine registration via WsiReg using shared morphology markers

CellScape captures quantitative FOXP3 expression validated by same-tissue HCR™ Gold RNA-FISH




Fig. 3 Per-cell FOXP3 protein expression correlates strongly with HCR™ Gold RNA-FISH transcript abundance ($r=0.96$), validating quantitative CellScape protein measurements.

Quantitative Treg Stratification into FOXP3_{low}, FOXP3_{mid}, and FOXP3_{high} phenotyping




Fig. 4 Co-expression of nuclear FOXP3 and cytoplasmic CD4 CellScape protein expression at single-cell resolution stratifies Treg populations into three quantitative phenotypic tiers.

TGFβ3 transcript abundance decreases progressively from FOXP3_{low} to FOXP3_{high} Tregs




Fig. 5 FOXP3_{low} Tregs show the highest proportion of TGFβ3⁺ cells, decreasing progressively through FOXP3_{mid} to FOXP3_{high}.

TGFβ3 drives pro-inflammatory Th17 differentiation via IL-23/IL-6 signaling, linking FOXP3_{low} Tregs to immune instability rather than canonical suppression.

References

- Christians et. Al. Development of a novel tri-omic spatial biology assay that provides in-depth immune profiling of the tumor immune microenvironment. AACR (2025).
- Li et. al. Spatial transcriptomics atlas of inflammatory bowel disease to guide implementation in research consortiums and clinical trials. Nat. Comms (*in press*)

Conclusion

- Quantitative FOXP3 reveals that FOXP3_{low}, FOXP3_{mid}, and FOXP3_{high} Tregs occupy spatially and functionally distinct niches within the inflamed colon.
- FOXP3_{low} and FOXP3_{mid} Tregs, enriched in stromal-vascular, cytotoxic, and remodeling niches and biased toward TGFβ3, suggest a destabilized pro-inflammatory Treg state.
- Cross-platform spatial multi-omics uncovers context-dependent Treg heterogeneity invisible to either modality alone.

CellScape™ PSP and CosMx® SMI

Scan here to
download or
learn more

