

# Linking single-cell spatial multiomics of TCR repertoires using CosMx SMI with live T cell function using the Beacon to characterize immune cell dynamics in tumors

Margaret L. Hoang<sup>1</sup>, Rachel Liu<sup>1</sup>, Dan McGuire<sup>1</sup>, Ashley Heck<sup>1</sup>, Kim Young<sup>1</sup>, Megan Vanderberg<sup>1</sup>, Claire Williams<sup>1</sup>, Martin Shelton<sup>1</sup>, Prajan Divakar<sup>1</sup>, Patrick Danaher<sup>1</sup>, Erin Piazza<sup>1</sup>, Joseph Beechem<sup>1</sup>

<sup>1</sup>Bruker Spatial Biology, Seattle, Washington, USA

## Introduction

Understanding how immune cells interact with tumor cells and the surrounding microenvironment is critical for predicting immunotherapy response. High-resolution spatial mapping of tumor-infiltrating lymphocytes (TILs) enables direct visualization of these interactions within tumor architecture. Using the CosMx<sup>®</sup> Spatial Molecular Imaging (SMI) platform, we integrate T-cell receptor (TCR) repertoire profiling with Whole Transcriptome (WTX) and immuno-oncology protein readouts to characterize immune phenotypes in FFPE tumor tissues. This approach captures TIL clonotypes and functional states, including activation, exhaustion, and exclusion, in situ, providing a detailed view of immune niches and their relationship to tumor and stromal features. Together, data highlight spatially resolved immune dynamics that inform biomarker discovery and therapeutic response.

## TCR diversity and CosMx TCR detection




Fig. 1 Each T-cell expresses a single, unique TCR generated through genomic recombination during development (a). TCRs are made by cutting and pasting individual segments together- V (Variable), J (Joining), and C (Constant) are joined for alpha chains; V, D (Diversity), J, and C for beta chains (a). The CosMx TCR assay detects these segments (panel content in c) alongside the whole transcriptome using barcoded probes and several rounds of reporter binding and imaging (b). Cells are segmented using protein markers and decoded RNA signals are assigned at single-cell resolution.

## CosMx TCR assay is specific in cell lines and FFPE tissue






Fig. 2. a) CosMx TCR panel measurements show strong concordance with RNAseq data in the CCRF-CEM T-cell line and detect the expected enrichment of the TRBV3-1 segment, supporting the specificity of the TCR assay. B) TCR probe signal is enriched in T-cells relative to other cell types and background across head and neck cancer FFPE tissue, confirming the specificity of TCR panel targets in FFPE tissue context.

## Linking live T-cell function to spatial context in head and neck cancer via CosMx multiomics and TCR repertoire profiling



## CosMx WTX + TCR + protein finds tumor-killing T-cells in space, exposing mechanisms of immune evasion by tumors



## Conclusion

- CosMx SMI enables single-cell, spatial TCR segment and clonotype identification with integrated WTX and protein in FFPE tumors.
- Spatial TCR mapping identifies T-cells shown to functionally kill tumors in their spatial context
- Combined WTX, high-plex protein, and TCR profiling offers a comprehensive view of the tumor microenvironment surrounding functional T-cells, exposing mechanisms of immune evasion and guiding patient specific therapeutic strategies

References:  
Zenga, J., Awan, M.J., Frei, A. et al. Radiation therapy results in preferential tumor antigen-specific lymphodepletion in head and neck cancer. *Nat Commun* 16, 5660 (2025).  
<https://doi.org/10.1038/s41467-025-60827-w>

Scan here to download or learn more

