

# SPACE: Spatially resolved multiomic analysis for high-throughput CRISPR screening in 3D models

Yi Cui<sup>1</sup>, Mengwei Hu<sup>2</sup>, Qianhui Huang<sup>2</sup>, Khoi Chu<sup>2</sup>, Sierra McKinzie<sup>1</sup>, Michael Patrick<sup>1</sup>, Sharanya Iyengar<sup>2</sup>, Maerjiangan Abuduli<sup>2</sup>, Marianne Spatz<sup>2</sup>, Nandita Joshi<sup>2</sup>, Brendan Miller<sup>2</sup>, Shams Vellarikkal<sup>2</sup>, Timothy Riordan<sup>1</sup>, Danny Bitton<sup>2</sup>, Jan Lubojacky<sup>2</sup>, Iya Khalil<sup>2</sup>, Federica Piccioni<sup>2</sup>, Michael Rhodes<sup>1</sup>, Alex Tamburino<sup>2</sup>, Shanshan He<sup>1</sup>, Joseph Beechem<sup>1</sup>, Vanessa Peterson<sup>2</sup>


1. Bruker Spatial Biology, Seattle, WA, USA; 2. Merck & Co, Inc., Rahway, NJ, USA




## Introduction

High-content single-cell perturbation screens are pivotal for elucidating gene functions and uncovering novel biology, yet conventional methods necessitate cell dissociation, forfeiting critical spatial information essential for dissecting cell-cell interactions and tissue architecture in complex microenvironments. While spatial CRISPR screening mitigates this partially, existing technologies are constrained by hypothesis-driven phenotyping panels limited to sparse RNA or protein coverage, curtailing comprehensive gene function assessment and discovery breadth. To overcome these barriers, we developed **SPACE (SPatial Cell Exploration)** assay, a pioneering platform that fuses whole-transcriptome profiling, CRISPR perturbations (sgRNAs and associated Unique Guide Identifiers UGIs), and multiplexed protein detection at single-cell resolution within intact 3D tissue contexts. SPACE delivers unbiased, transcriptome-wide readouts alongside 68 protein markers (64 CosMx Human IO protein panel targets plus 4 morphological markers), vastly expanding phenotypic landscapes in spatial screens. As the highest-plex coverage multimodal CRISPR assay to date, SPACE achieves this at unprecedented scale and affordability – profiling hundreds of spheroids encompassing tens of thousands of cells per slide – largely outperforming sequencing-based alternatives in efficiency.

## SPACE Assay Design




**Fig. 1** Diagram of SPACE experimental workflow. Primary cancer associated fibroblasts (CAF) are CRISPR edited, co-cultured with BxPC3 cancer cells in ultra-low attachment plates to allow spheroid formation, and fixed. The spheroids are then embedded in FFPE blocks while maintaining 96-well arrayed format for CRISPR identity confirmation. The sectioned slides undergo SPACE profiling of CRISPR molecules, transcriptome and protein on the same slide




**Fig. 2** Example images from SPACE workflow. Left: morphological staining of whole slide containing a spheroid array. Top right: Morphological staining markers. Bottom right: composite images of representative spheroids, showing different synthetic sgRNAs, COL1A1 mRNA, PanCK protein and nuclear staining.

## SPACE Assay Characterization




**Fig. 3** SPACE is compatible with different panel sizes, showing consistent detection of CRISPR targets with high specificity. Left: Summary table for SPACE performance in RNA detection. Middle: Correlation of detected CRISPR molecules (upper: CRISPR sgRNAs; lower: UGIs) count between WTX+CRISPR (x axis) and 6K+CRISPR (y axis). Right: Counts per CAF cell of the top 5 most abundant sgRNAs and UGIs in AAVS1-KO and ISG20-KO spheroids.




**Fig. 4** Primary cell identification. Left: cell clustering using RNA data. Middle: Differentially expressed genes across main clusters. Right: Cell annotation with key marker genes.

## SPACE Revealed Novel Biomarkers in CAF




**Fig. 5** In ISG20-KO CAF, SPACE identified several MMPs significantly downregulated, which were later validated with orthogonal assays at both RNA and protein levels.




**Fig. 6** Kaplan-Meier plots of patient overall survival in different tumor types categorized as high-ISG20 or low-ISG20 subgroups. (Data source: TCGA/GTEX)

## SPACE Enabled Spatially Informed LR Analysis




**Fig. 7** Left: Score differences between the spatially informed and nonspatially informed LR analysis for each CRISPR edits (blue – high LR interaction score in nonspatially informed analysis; red – high LR interaction score in spatially informed analysis). Only LR pairs identified in more than 5 KO conditions are presented. Right: Chord diagrams showing the spatially informed LR pairs specific to ISG20-KO and RNF213-KO spheroids compared to nonspatially informed analysis (LR pairs with score difference > 0).


## Same-Sample Multimodal SPACE



**Fig. 8** Diagram of multimodal SPACE workflow. Samples were first stained with CosMx 64-plex protein panel, followed by the SPACE WTX workflow for RNA profiling.




**Fig. 9** Representative images of protein detected in multimodal SPACE.




**Fig. 10** Left-Middle: Multiomic SPACE revealed the same cell clustering outcomes as WTX-only SPACE (Fig. 3). Right: Representative image of RNA, protein and CRISPR targets detected in the same cells.

## SPACE for Spatial Phenotyping



**Fig. 11** Tumor cell EMT score calculated using protein data per CRISPR editing condition



**Fig. 12** Left: Spatially resolved heatmap of EMT score in each individual tumor cells, with 4 spheroids zoomed in (Upper right). Lower right: Tumor cell EMT score negatively correlates with their distance to the closest CAF cell.

## Summary

- SPACE integrates whole transcriptome, multiplexed proteins, and CRISPR sgRNA ID at subcellular resolution in 3D models, providing unparalleled phenotypic depth.
- SPACE scales to hundreds of spheroids per slide for cost-effective screening.
- SPACE is FFPE-compatible and supports multimodal, orthogonal cross-validation.
- SPACE preserves spatial context to reveal new regulators and mechanisms.
- SPACE accelerates drug discovery with disease-relevant functional genomics datasets.

**Reference:**  
 bioRxiv 2025.09.14.675819;  
 doi: <https://doi.org/10.1101/2025.09.14.675819>  
 bioRxiv 2024.11.27.625536;  
 doi: <https://doi.org/10.1101/2024.11.27.625536>

Scan here to  
 download or  
 learn more

