

## ChipCytometry™ ZellSafe™ FFPE Tissue Chip Kits

ZellSafe™ FFPE Tissue Chip Kits contain chips and reagents needed to load and store samples for ChipCytometry™ assays using FFPE tissues and the ZellScannerONE™ Instrument. Store loaded tissue chips at 2-8°C and unloaded chips at room temperature.

ZellSafe™ Chip Kits for FFPE Tissues			
Item	Size	Catalog #	Contents
ZellSafe™ FFPE Tissue Chip Kit	10 chips	28050606/04-010	10 ZellSafe™ Chips 40 Pre-treated Coverslips 1 x 60 mL ZKW Wash Buffer 1 x 50 mL ZKW Storage Buffer 1 Chip Storage Box
ZellSafe™ FFPE Tissue Chip Kit	50 chips	28050606/04-050	50 ZellSafe™ Chips 200 Pre-treated Coverslips 1 x 500 mL ZKW Wash Buffer 3 x 50 mL ZKW Storage Buffer 2 Chip Storage Boxes
ZellSafe™ FFPE Tissue Chip Kit	100 chips	28050606/04-100	100 ZellSafe™ Chips 400 Pre-treated Coverslips 2 x 500 mL ZKW Wash Buffer 6 x 50 mL ZKW Storage Buffer 3 Chip Storage Boxes

To learn more, visit [canopybiosciences.com/chipcytometry](https://canopybiosciences.com/chipcytometry) or email us at [info.canopy@bruker.com](mailto:info.canopy@bruker.com)

**For Research Use Only.** Not for use in diagnostic procedures.

© 2022 Canopy Biosciences. All rights reserved. All trademarks are the property of Canopy Biosciences, a Bruker Company unless otherwise specified.

## ChipCytometry™ ZellSafe™ Fresh Frozen Tissue Chip Kits

ZellSafe™ Fresh Frozen Tissue Chip Kits contain chips and reagents needed to load and store samples for ChipCytometry™ assays using fresh frozen tissues and the ZellScannerONE™ Instrument. Store loaded tissue chips at 2-8°C and unloaded chips at room temperature.

ZellSafe™ Chip Kits for Fresh Frozen Tissues			
Item	Size	Catalog #	Contents
ZellSafe™ Fresh Frozen Tissue Chip Kit	10 chips	28050606/02-010	10 ZellSafe™ Chips 10 Coverslips 1 x 60 mL ZKW Wash Buffer 1 x 50 mL ZKW Storage Buffer 1 x 60 mL ZKW Fixation Buffer 1 Chip Storage Box
ZellSafe™ Fresh Frozen Tissue Chip Kit	50 chips	28050606/02-050	50 ZellSafe™ Chips 50 Coverslips 1 x 500 mL ZKW Wash Buffer 3 x 50 mL ZKW Storage Buffer 3 x 60 mL ZKW Fixation Buffer 2 Chip Storage Boxes
ZellSafe™ Fresh Frozen Tissue Chip Kit	100 chips	28050606/02-100	100 ZellSafe™ Chips 100 Coverslips 2 x 500 mL ZKW Wash Buffer 6 x 50 mL ZKW Storage Buffer 6 x 60 mL ZKW Fixation Buffer 3 Chip Storage Boxes

To learn more, visit [canopybiosciences.com/chipcytometry](https://canopybiosciences.com/chipcytometry) or email us at [info.canopy@bruker.com](mailto:info.canopy@bruker.com)

**For Research Use Only.** Not for use in diagnostic procedures.

© 2022 Canopy Biosciences. All rights reserved. All trademarks are the property of Canopy Biosciences, a Bruker Company unless otherwise specified.

## ChipCytometry™ ZellSafe™ Cell Chip Kits

ZellSafe™ Cell Chip Kits contain chips and reagents needed to load and store samples for ChipCytometry™ assays using cell suspensions and the ZellScannerONE™ Instrument. Store chips for cell suspensions at 2-8°C.

ZellSafe™ Chip Kits for Cell Suspensions			
Item	Size	Catalog #	Contents
ZellSafe™ Cell Chip Kit	10 chips	28050606/01-010	10 ZellSafe™ Chips 1 x 60 mL ZKW Wash Buffer 1 x 50 mL ZKW Storage Buffer 1 x 60 mL ZKW Fixation Buffer 1 Chip Storage Box
ZellSafe™ Cell Chip Kit	50 chips	28050606/01-050	50 ZellSafe™ Chips 3 x 60 mL ZKW Wash Buffer 1 x 50 mL ZKW Storage Buffer 1 x 60 mL ZKW Fixation Buffer 2 Chip Storage Boxes
ZellSafe™ Cell Chip Kit	100 chips	28050606/01-100	100 ZellSafe™ Chips 1 x 500 mL ZKW Wash Buffer 3 x 50 mL ZKW Storage Buffer 3 x 60 mL ZKW Fixation Buffer 3 Chip Storage Boxes

To learn more, visit [canopybiosciences.com/chipcytometry](https://canopybiosciences.com/chipcytometry) or email us at [info.canopy@bruker.com](mailto:info.canopy@bruker.com)

**For Research Use Only.** Not for use in diagnostic procedures.

© 2022 Canopy Biosciences. All rights reserved. All trademarks are the property of Canopy Biosciences, a Bruker Company unless otherwise specified.

## ChipCytometry™ ZellSafe™ Rare Cell Chip Kits

ZellSafe™ Rare Cell Chip Kits contain chips and reagents needed to load and store samples for ChipCytometry™ assays using cell suspensions and the ZellScannerONE™ Instrument. Store chips for cell suspensions at 2-8°C.

ZellSafe™ Chip Kits for Rare Cells			
Item	Size	Catalog #	Contents
ZellSafe™ Rare Cell Chip Kit	10 chips	28050606/03-010	10 ZellSafe™ Chips 1 x 60 mL ZKW Wash Buffer 1 x 50 mL ZKW Storage Buffer 1 x 60 mL ZKW Fixation Buffer 1 Chip Storage Box
ZellSafe™ Rare Cell Chip Kit	50 chips	28050606/03-050	50 ZellSafe™ Chips 1 x 500 mL ZKW Wash Buffer 3 x 50 mL ZKW Storage Buffer 3 x 60 mL ZKW Fixation Buffer 2 Chip Storage Boxes
ZellSafe™ Rare Cell Chip Kit	100 chips	28050606/03-100	100 ZellSafe™ Chips 2 x 500 mL ZKW Wash Buffer 6 x 50 mL ZKW Storage Buffer 6 x 60 mL ZKW Fixation Buffer 3 Chip Storage Boxes

To learn more, visit [canopybiosciences.com/chipcytometry](https://canopybiosciences.com/chipcytometry) or email us at [info.canopy@bruker.com](mailto:info.canopy@bruker.com)

**For Research Use Only.** Not for use in diagnostic procedures.

© 2022 Canopy Biosciences. All rights reserved. All trademarks are the property of Canopy Biosciences, a Bruker Company unless otherwise specified.

NOV 2022 V2